

Perfetta®

Children Face Mask



PHU BAO GROUP

**HIGH
QUALITY**



- ✓ For Children Use (6-12 yrs)
- ✓ Procedural mask
- ✓ Moderate Barrier Protection
- ✓ BFE at 0.3 micron >99%
- ✓ PFE at 0.1micron >99%
- ✓ Soft & Comfortable
- ✓ Breathable suitable for Children

PRODUCT DESCRIPTION

Product Name	<ul style="list-style-type: none"> ▪ Perfetta Children Face Mask, 3Ply Single-use.
Model Number	<ul style="list-style-type: none"> ▪ PERTR3W50 (White color) PER3B50 (Blue color)
US Standard	<ul style="list-style-type: none"> ▪ Masks tested at Nelson Lab, USA according to ASTM F2100-19 ▪ BFE ≥ 99%, Delta P <5 mmH2O/cm2, PFE ≥ 99%,
EU Standard	<ul style="list-style-type: none"> ▪ Masks tested at Nelson Lab, USA according to EN14683:2019
Intended Use	<ul style="list-style-type: none"> ▪ Perfetta Children– 3Ply Single-use Face Mask covers the user’s nose and mouth and meet the certain particulate filtration efficiency to prevent patient exposure to large particles. ▪ This is a single-use, disposable, provided non-sterile. Pleat O3 fold, Ear loop type mask.
Mask Dimension	<ul style="list-style-type: none"> ▪ 145mm x 95 mm (+-5mm)
Materials	<ul style="list-style-type: none"> ▪ 3 Layers composite of: <ul style="list-style-type: none"> - Outer Layer: Polypropylene spunbond non-woven fabric. - Middle Layer: Polypropylene melt blown filter nonwoven fabric. - Inner Layer: Soft Polypropylene spunbond non-woven fabric. ▪ Nose clip: 100% PE ,Size 80 mm (+- 5mm). Color White. ▪ Ear loop: polyester-spandex elastic. Latex-Free. Width 3.5mm, Length 150 mm. White color

PACKAGE

Expiration date	<ul style="list-style-type: none"> ▪ 3 years from production date. Please see the manufacturing date in the box (MFG date)
Packaging	<ul style="list-style-type: none"> ▪ 35 pcs Individually Wrapped /box, 40 boxes/carton, 1,400 pcs/carton
Carton / Case	<ul style="list-style-type: none"> ▪ Size: 520mmx 330mmx 420mm, vertical packing.
Storage	<ul style="list-style-type: none"> ▪ Cool & dry place, avoid direct sunlight and high temperatures. The humidity of no more than 80%, away from corrosive gas in a well-ventilated area. Avoid high temperatures.
Warning	<ul style="list-style-type: none"> ▪ This product is not to be used for anti-microbial protection, antiviral protection and it does not eliminate or decrease the risk of contracting any disease of infection. ▪ One time use only. Please dispose mask after use. ▪ Do not use if the packaging is damaged. ▪ Please stop wearing if get sensitive or allergic. ▪ This product is distributed under FDA’s Enforcement Policy for Face Masks & Respirators during Coronavirus Disease (COVID-19) Public Health Emergency.

Remark: All information presenting in this document is base on the latest data available at the time of printing, Phu Bao Medical reserve the right to make change at anytime without prior notice.