Single-Use Medical Face Mask



PHU BAO GROUP





 Perfetta Moderate – 3Ply Single-use Medical Face Mask 	
LADIESPER3W50 (White color)	
 Masks tested at Nelson Lab, USA according to ASTM F2100-19 Level 2 BFE ≥ 99%, PFE ≥ 99%, Delta P <6 mmH20/cm2, Fluid resistance 120mm Hg, Flammability Class 1. Material meet Biocompatible ISO10993 : Non-Irritation, Non-Sensitization, Non-Cytotoxicity. 	
 Masks tested at Nelson Lab, USA according to EN14683:2019 Type IIR. 	
 Perfetta Moderate – 3 Ply Single-use Medical Face Mask covers the user's nose and mouth and meet the certain fluid barrier, flammability and particulate filtration efficiency to prevent Healthcare Personal or patient exposure to respiratory droplets and large particles. This is a single-use, disposable, provided non-sterile. Pleat 03 fold, Ear loop type mask. 	
 150 mm x 95 mm (+-5mm) 	
 3 Layers composite of: Outer Layer: Polypropylene spunbond non-woven fabric. Middle Layer: Polypropylene melt blown filter nonwoven fabric. Inner Layer: Soft Polypropylene spunbond non-woven fabric. Nose clip: 100% PE OR metal wire with plastic covering, Size 90mm (+- 5mm). Color White. Ear loop: polyester-spandex elastic. Latex-Free. Width 3.5mm, Length 170mm. White color 	
PACKAGE	
 3 years from production date. Please see the manufacturing date in the box (MFG date) 	
 50 pcs/box, 50 boxes/carton, 2,500 pcs/carton 	
Size: 535mmx 340mmx 470mm, vertical packing.	
 Cool & dry place, avoid direct sunlight and high temperatures. The humidity of no more than 80%, away from corrosive gas in a well-ventilated area. Avoid high temperatures. 	
 This product is not to be used for anti-microbial protection, antiviral protection and it does not eliminate or decease the risk of contracting any disease of infection. One time use only. Please dispose mask after use. Do not use if the packaging is damaged. Please stope wearing if get sensitive or allergic. This product is distributed under FDA's Enforcement Policy for Face Masks & Respirators during Coronavirus Disease (COVID-19) Public Health Emergency. 	

Remark: All information presenting in this document is base on the latest data available at the time of printing, Phu Bao Medical reserve the right to make change at anytime without prior notice.