

Perfetta Ultra®

Single-Use Medical Face Mask



PHU BAO GROUP

LEVEL 3
ASTM



- ✓ Medical & Professional Use
- ✓ High Barrier Protection
- ✓ Fluid Resistant 160mm Hg
- ✓ PFE at 0.1micron >99%
- ✓ BFE at 0.3 micron >99%
- ✓ VFE at 3 micron > 99%
- ✓ Hypoallergenic

PRODUCT DESCRIPTION

Product Name	<ul style="list-style-type: none"> ▪ Perfetta Ultra - 4Ply Single-use Medical Face Mask
Model Number	<ul style="list-style-type: none"> ▪ ULT4W50 (White color) ULT4B50 (Blue color)
US Standard	<ul style="list-style-type: none"> ▪ Masks tested at Nelson Lab, USA according to ASTM F2100-19 Level 3 ▪ BFE, PFE, VFE ≥ 99%, Delta P <6 mmH2O/cm2, Fluid resistance 160mm Hg, Flammability Class 1. ▪ Meet Biocompatible ISO10993 : Non-Irritation, Non-Sensitization, Non-Cytotoxicity.
EU Standard	<ul style="list-style-type: none"> ▪ Masks tested at Nelson Lab, USA according to EN14683:2019 Type IIR.
Intended Use	<ul style="list-style-type: none"> ▪ Perfetta Ultra – 4 Ply Single-use Medical Face Mask covers the user’s nose and mouth and meet the certain fluid barrier, flammability and particulate filtration efficiency to prevent Healthcare Personal or patient exposure to respiratory droplets and large particles. ▪ This is a single-use, provided non-sterile. Pleat 04 fold Omega shape, Ear loop type mask.
Mask Dimension	<ul style="list-style-type: none"> ▪ 175mm x 95 mm (+-5mm)
Materials	<ul style="list-style-type: none"> ▪ 4 Layers composite of: <ul style="list-style-type: none"> - Outer Layer: Polypropylene spunbond non-woven fabric. Color white / Blue. - 1st Middle Layer: Polypropylene spunbond non-woven fabric. Color White. - 2nd Middle Layer: Polypropylene melt-blown filter nonwoven fabric. Color White. - Inner Layer: Soft Polypropylene spunbond non-woven fabric. Color White. ▪ Nose clip: 100% PE OR metal wire with plastic covering. Size 90mm (+ - 5mm). Color White. ▪ Ear loop: polyester-spandex elastic. Latex-Free. Width 3.5mm, Length 160mm. White color

PACKAGE

Expiration date	<ul style="list-style-type: none"> ▪ 3 years from production date. Please see the manufacturing date in the box (MFG date)
Packaging	<ul style="list-style-type: none"> ▪ 50 pcs/box, 40 boxes/carton, 2,000 pcs /carton
Carton / Case	<ul style="list-style-type: none"> ▪ Size: 555mmx 400mmx 470mm, vertical packing.
Storage	<ul style="list-style-type: none"> ▪ Cool & dry place, avoid direct sunlight and high temperatures. The humidity of no more than 80%, away from corrosive gas in a well-ventilated area. Avoid high temperatures.
Warning	<ul style="list-style-type: none"> ▪ This product is not to be used for anti-microbial protection, antiviral protection and it does not eliminate or decrease the risk of contracting any disease of infection. ▪ One time use only. Please dispose mask after use. ▪ Do not use if the packaging is damaged. ▪ Please stop wearing if get sensitive or allergic. ▪ This product is distributed under FDA’s Enforcement Policy for Face Masks & Respirators during Coronavirus Disease (COVID-19) Public Health Emergency.

Remark: All information presenting in this document is base on the latest data available at the time of printing, Phu Bao Medical reserve the right to make change at anytime without prior notice.